**Praxiject™ Normal Saline**  
 Prefilled Syringe with 0.9% Sodium Chloride Injection USP  
 Sterile – Non Pyrogenic – Single Use

**INDICATIONS AND CLINICAL USE**  
 For vascular access devices flushing only. May be placed on sterile field.  
**CONTRAINDICATIONS**  
 Not for dry product reconstitution or for medication dilution.

**WARNINGS AND PRECAUTIONS**  
 Use aseptic technique. Single use only. Do not reuse. Not for direct intravenous injection. Do not use if package is opened or damaged. Do not use if the cap/adaptor on the syringe is not intact. Visually inspect the contents of each prefilled syringe for clarity, particulate matter, precipitate, discoloration and leakage prior to use. Do not use if any of the aforementioned is observed. Do not allow air to be trapped in fluid path. Do not resterilise. Sterile unless package opened or damaged.

**ADVERSE REACTIONS**  
 No known adverse reactions when the product is used as indicated.

**INSTRUCTIONS FOR USE**  
**USE ASEPTIC TECHNIQUE.**

- 1 Aseptically open package containing the 0.9% Sodium Chloride prefilled syringe. Visually inspect product for particulate matter and discoloration prior to use according to WARNINGS AND PRECAUTIONS.
- 2 Remove protector as required. Expel air from syringe. Do not allow air to be trapped in fluid path. Use only with compatible luer lock connectors.
- 3 Remove and discard unused portions and empty syringes according to the institution's biohazardous waste policy.

**ACTION**  
 0.9% Sodium Chloride Injection USP solution is a sterile, aqueous solution having approximately the same osmotic pressure and composition as extracellular fluids. It is non-irritating to tissues. It is used to flush vascular access devices in order to maintain catheter patency and to prevent contact between incompatible medications or fluids.

**PRODUCT DESCRIPTION AND PACKAGING**  
 Each syringe is supplied as a clear, latex-free, plastic, prefilled syringe, packaged in a heat-sealed pouch. Each syringe in an unopened, undamaged package contains 0.9% Sodium Chloride Injection USP with no preservatives.

Catalogue Number	Description	Packaging
3704	5mL 0.9% Sodium Chloride Injection USP in a 5mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 120 pouches (120 syringes) per product case
37043	3mL 0.9% Sodium Chloride Injection USP in a 5mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 120 pouches (120 syringes) per product case
3705C	10mL 0.9% Sodium Chloride Injection USP in a 10mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 100 pouches (100 syringes) per product case
37053	3mL 0.9% Sodium Chloride Injection USP in a 10mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 100 pouches (100 syringes) per product case
37055	5mL 0.9% Sodium Chloride Injection USP in a 10mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 100 pouches (100 syringes) per product case
3706	20mL 0.9% Sodium Chloride Injection USP in a 20mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 60 pouches (60 syringes) per product case

**STORAGE AND STABILITY**  
 0.9% Sodium Chloride Solution Injection USP should be stored between 15°C and 30°C (59°F and 86°F), protected from direct sunlight and freezing. Product may be stored as indicated by the expiration date on the product package.

**SYMBOLS ON PRODUCT LABELS**

REF: Catalogue (Reference) Number	Caution, consult accompanying documents	STERILE R Sterile by irradiation
LOT Lot number	No latex	Do not reuse, single use only
Use by	Recommended Storage Temperature	Authorised representative in the European community Emergo-Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
CE marking	Fill volume	

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