

NOVA BETAFIX[®] AG

Elastic fistula bandage provided with non-adherent ionic silver based pad and fixing plaster



Description

Elastic fistula bandage with non-adherent silver ions based pad provided with plaster for better tightness.

Intended use

Compressive hemostatic action, specific for post-dialysis AVF management to be used after the removal of fistula needles. Also suitable for patients who require intravenous treatment and blood transfusions.

Indication

The presence of the silver-based layer in **NOVA BETAFIX[®] AG** carries an immediate and long-term antibacterial activity, helping to reduce the risk of infection associated with the use of intravascular devices.

Available Models

CODE	LENGHT	WIDTH	SILVER PAD	PCS/PACK	PACKS/BOX
S0014/AG	45 cm ±10%	6 cm	3.5 x 3.5 cm	2	50
S0015/AG				1	
S0029/AG			2.5 x 2.5 cm	1	100
S0030/AG				2	50

⚠ Do not use silver dressings on patients sensitive to silver

Made in Italy by
EMODIAL

SILVER LINE

