

BETAFIX®

Elastic fistula bandage provided with a two layers highly absorbent pad and plaster for a better tightness



Description

Elastic fistula bandage provided with a two layers highly absorbent pad and adhesive for a better tightness. Also suitable for patients who require intravenous treatment and blood transfusion.

Intended use

Compressive haemostatic action, specific for post dialysis arteriovenous access care. To be used after the removal of the fistula needles. Also suitable for patients who require intravenous treatment and blood transfusion.

Indication

Designed to stop bleeding and to reduce the hemorrhagic risk, thanks to the compression on the needle holes, also working as a barrier against the risk of infections. The non-adherent, biocompatible layer on the pad doesn't stick to the skin avoiding the clot removal and the risk of bleeding. The cotton bandage avoids the sensitization due to the glue of adhesive.

Available Models

CODE	LENGHT	WIDTH	PAD	PCS/PACK	PACKS/BOX
S0029	50 cm ± 10%	6 cm	ø 2.6 cm	1	500
S0029-1				2	
S0030			ø 3.5 cm	300	
S0041			3.5 x 4.5 cm		
S0014	55 cm ± 10%			1	400
S0029/A				ø 2.6 cm	
S0030/A	25 cm ± 10%			2	300
S0030/R					

Made in Italy by
 **EMODIAL**

STANDARD LINE