


EXIT-PAD[®] AG

Exit-site antiseptic silver dressings



Description

Non adherent silver ions based dressing designed to allow the catheter positioning.

Intended use

Dressing of the points of insertion of intravascular devices. To be used as a mechanical barrier for compression and absorption of exudates Compatible with all types of catheter use.

Indication


The presence of the silver-based layer carries an immediate and long term antibacterial action. It helps to reduce the risk of local infection of the exit-site, of the tunnel or systemic (catheter-correlated sepsis) related to the use of intravascular devices. It reduces the signs of inflammation and promotes tissue repair.

Available Models


CODE	MODEL	HOLE	SIZE	PCS/PACK	PACKS/BOX
S0903/AG	Cut and Hole	ø 4 mm	3.5 x 4.5 cm	1	100
S0904/AG	Cut and Hole	ø 5 mm	3.5 x 4.5 cm	1	100
S0905/AG	Y-Cut		5 x 4 cm	1	100
S0925/AG	Cut and hole	ø 4 mm	3.5 x 2.2 cm	2	100

⚠ Do not use silver dressings on patients sensitive to silver

Made in Italy by
EMODIAL



CUT AND HOLE



Y CUT *

* Y CUT model is suitable for PICC line catheters

SILVER LINE

