

EXIT-POCKET® AG

Exit-site protective pocket with antiseptic silver pad



Description

Non-woven pocket closed on three sides provided with a silver ion-releasing pad. Available on two supports: on breathable hypoallergenic non-woven (NW) and on waterproof polyurethane membrane (PU)

Intended use

The peculiar configuration of the device allows two functions:

- active dressing of the exit-site thanks to the silver ion-releasing pad
- containment of central venous or peritoneal catheter, to protect the skin around the exit-site.

Indication

The presence of the silver-based layer carries an immediate and long term antibacterial action. It helps to reduce the risk of local infection of the exit-site, of the tunnel or systemic (catheter-correlated sepsis) related to the use of intravascular devices.

Reduces the signs of inflammation and promotes tissue regeneration. The containment pocket also prevents redness of the skin caused by the mechanical action of rubbing the catheter.

Available Models

CODE	POCKET	SILVER PAD	SUPPORT	CATHETER	VERSION	PCS/PACK	PACKS/BOX	
S2021YH/AG	21 x 11 x 8 cm	5 x 4 cm	PU	TESIO	Opening/ Closing System	1	50	
S2021Y/AG			NW					
S2027YH/AG	PU							
S2027Y/AG	NW							
S2021FH/AG	21 x 13 x 8 cm		PU	MONO LUMEN				
S2021F/AG			NW					
S2027FH/AG	27 x 19 x 8 cm		PU					MONO LUMEN
S2027F/AG			NW					
S2024FH/AG	24 x 18 x 8 cm		PU	MONO LUMEN	Open			
S2030FH/AG	30 x 24 x 8 cm							

⚠ Do not use silver dressings on patients sensitive to silver

Made in Italy by
EMODIAL



SILVER LINE