

# EXIT-POCKET® PT

Nonwoven pocket  
provided with  
elastic belt  
and Velcro  
fastening/  
closing system



#### Description

Nonwoven pocket provided with elastic belt and Velcro fastening/closing system.

#### Intended use

Contains the catheter lumen between each dialysis treatments.

#### Indication

Protects catheter surrounding skin and acts as a barrier against microbial contaminations.

#### Available Models

CODE	POCKET SIZE	PCS/PACK	PACKS/BOX
S2013/PT	18 x 8 cm	1	40
S2010/PT	24 x 8 cm		

Made in Italy by  
**EMODIAL**

STANDARD LINE