

# EXIT-POCKET®

Exit-site  
dressing device  
with non adherent  
pad set on a  
protective pocket



## Description

Non-Woven pocket closed on three sides for CVC protection. Available with or without Y-Cut for catheter exit.

## Intended use

Protection of the exit-site surrounding skin by containing the lumen(s) of the catheter in between dialysis sessions. Non-invasive device intended to use as a mechanics barrier for the compression and absorption of exudates.

## Indication

Containment of bleeding and exudates. Prevents skin redness and microbial contamination of the exit-site.

## Available Models

| CODE       | POCKET         | VERSION | PCS/PACK | PACKS/BOX |
|------------|----------------|---------|----------|-----------|
| S2010A/P1P | 30 x 24 x 8 cm | Y-Cut   | 50       | 300       |
| S2011A/P1P | 18 x 12 x 8 cm |         |          |           |
| S2013A/P1P | 24 x 18 x 8 cm |         |          |           |
| S2010Y/P1P | 30 x 24 x 8 cm |         |          |           |
| S2011Y/P1P | 18 x 12 x 8 cm |         |          |           |
| S2013Y/P1P | 24 x 18 x 8 cm |         |          |           |

Made in Italy by  
**EMODIAL**

**STANDARD LINE**