




Fistula plaster  
provided with  
antiseptic  
silver pad on  
hypoallergenic  
adhesive



### Description

Not adherent pad with silver ion releasing layer applied on breathable non-woven fabric patch coated with a silicone-based hypoallergenic adhesive gel with **SoftTouch** technology.

### Intended use

Compressive hemostasis in extracorporeal and infusion treatments. Particularly suitable for the management of patients who need frequent dressings of the needles insertion points and also are subject to MARS (Medical Adhesive Related Skin Injuries) problems, such as dialysis patients.

### Indication

The dressing performs multiple actions:

- thanks to the presence of the silver ion-releasing film which has a broad-spectrum antimicrobial action, it helps to reduce the risk of infection of the venipuncture areas such as the arteriovenous fistula (AVF), the points of injection for blood samples and for the placement of peripheral venous accesses aimed at the administration of drugs or blood derivatives;
- reduces the signs of inflammation related to continuous venipunctures;
- helps to prevent the formation of eschar, making subsequent venipuncture less traumatic for the patient;
- the **SoftTouch** adhesive film preserves the skin from mechanical damage caused by the continuous removal of the patch, a consequence of continuous treatments, favoring tissue repair.

Optimal solution for patients who experience skin conditions with dry, moist or scaly skin, widespread conditions in patients undergoing, for example, dialysis treatment.

### Available Models

CODE	NW PLASTER	SILVER PAD	PCS/PACK	PACKS/BOX
S0074/AG	5 x 7.2 cm	2.8 x 2.8 cm	1	100

⚠ Do not use silver dressings on patients sensitive to silver

