


NOVA SOVAN[®] AG

Antiseptic breathable adhesive dressing provided with non-adherent ionic silver based pad



Description

Fistula pad provided with silver ion-releasing anti-adherent film. Set on breathable hypoallergenic non-woven plaster (NW) or waterproof polyurethane membrane (PU).

Intended use

Compressive hemostatic action, specific for post-dialysis arteriovenous access care to be used after the removal of fistula needles. Particularly suitable for the dressing of arteriovenous fistulas managed with the method of button-hole.

Indication

- Reduces the risk of the arteriovenous fistula infection thanks to the presence of the silver based film with long-term and immediate antibacterial action.
- Reduces signs of inflammation associated with repeated venipuncture of the AVF.
- Prevents eschar formation and makes less traumatic for the patient the next venipuncture.
- The waterproof adhesive model protects the fistula from contact with water during patient's hygiene.

Available Models

CODE	NW PLASTER	PU FILM	SILVER PAD	PCS/PACK	PACKS/BOX
S0075/AG	15 x 5 cm			2	100
S0079/AG	10 x 5 cm		2.5 x 2.5 cm	2	30
S0079-1/AG				1	40
S0079-1H/AG		10 x 5 cm		1	100

⚠ Do not use silver dressings on patients sensitive to silver

Made in Italy by
EMODIAL

SILVER LINE

