

SOVAN[®]

Highly absorbing
medical pad set
 on **hypoallergenic**
plaster designed
 to plug and stop
 bleeding at the
 end of the dialytic
 treatment



Description

Highly absorbing two layers hemostatic pad set on hypoallergenic, air permeable adhesive. Also suitable for patients who require intravenous treatment and blood transfusion.

Intended use

Compressive hemostatic action, specific for post dialysis arteriovenous access care. To be used for the hemostasis after the removal of the fistula needles.

Indication

Designed to stop bleeding and to reduce the hemorrhagic risk, thanks to the compression on the needle holes, also working as a barrier against the risk of infections. The non-adherent, biocompatible layer on the pad does not stick to the skin, avoiding the clot removal and the risk of re-bleeding. The adhesive helps to keep the pad in position.

Available Models

CODE	NW PLASTER	PAD	PCS/PACK	PACKS/BOX
S0074	15 x 5 cm	ø 2.6 cm	1	500
S0075			2	
S0079	10 x 5 cm			

Made in Italy by
 **EMODIAL**

STANDARD LINE